
ENGLISH P 5-8

ESPAÑOL P 9-12

# Celsite®

## Implantable Access Port System Sistema Implantable de Acceso

### Instructions for Use Instrucciones de uso

en (English): INSTRUCTIONS FOR USE  
**Celsite® IMPLANTABLE ACCESS PORT SYSTEM**  
 PLEASE READ CAREFULLY  
 TABLE OF CONTENTS

I - DESCRIPTION OF THE DEVICE	Page 2
II - PACKAGE CONTENTS	Page 4
III - INDICATIONS	Page 5
IV - CONTRAINDICATIONS	Page 5
V - POTENTIAL COMPLICATIONS (Immediate and delayed)	Page 5
VI - GENERAL WARNINGS & PRECAUTIONS	Page 5
VI-1 HANDLING AND STORAGE	
VI-2 PATIENT	
VI-3 CATHETER POSITION	
VI-4 DRUG INCOMPATIBILITIES	
VI-5 NEEDLES	
VI-6 IMPLANTATION TECHNIQUES	
VI-7 OCCLUSION	
VI-8 INFECTION	
VII - CENTRAL VENOUS IMPLANTATION TECHNIQUES	Page 6
VII-1 Venous catheters	
VII-2 Percutaneous Technique	
VII-3 Surgical Cut-down Technique	
VII-4 References requiring modifications to the above techniques	
VII-4-1 Brachial Implantation Technique	
VII-4-2 Braunule Technique	
VII-4-3 Valved Catheter	
VII-4-4 Pre-connected catheter and port	
VIII - USE OF THE PORT AND CATHETER	Page 7
VENOUS PORTS	
IX - MAINTENANCE OF THE PORT AND CATHETER	Page 7
IX-1 RINSING AND HEPARINIZATION	
VENOUS PORTS	
IX-2 NEEDLE WITHDRAWAL	
X - DURATION OF IMPLANTATION AND REMOVAL OF THE PORT	Page 8
X-1 DURATION OF IMPLANTATION OF THE SYSTEM	
X-2 REMOVAL OF THE SYSTEM	
X-3 SPECIAL ATTENTION: BABYPOR® S AND SMALL PORTS	
PRODUCT LIST	Page 13
FIGURES	Page 14
SYMBOLS GLOSSARY	Page 15

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en: The Celsite® Access Port Systems are implantable port and catheter systems that allow safe, repeated access to the patient's bloodstream. The port chamber and catheter design (See Figure A) can be used for the administration of medications and fluids. The Celsite® port systems consist of an access port with a silicone septum, which is connected to a catheter using a connection ring, and implantation accessories. The triangular shaped port has a low-profile nose, finger stops on the side of the housing, and a round base. Celsite® Access Ports have suture holes or suture zones to secure placement during implantation. The system is provided sterile and intended for single use only. Celsite® ports are offered in "standard" or "small" sizes. (A "standard" size has 0.50 or 0.25 mL dead volume of port; a "small" size has 0.15 mL dead volume of port.) The Celsite® access port family includes the below product lines. Refer to page 13 of this IFU for a complete list of Celsite® reference numbers and information regarding port/catheter materials, volume, implantation route, and implantation techniques.

- Celsite® Access Ports (Venous) are standard sized ports manufactured with a low-profile design with titanium/polysulphone/silicone or titanium/epoxy/silicone materials. Catheters are available in either polyurethane or silicone materials (REF: 05430425, 05432460, 05433823, 05433750, 05433556, 05436946, 05430095). Some Celsite® ports are available with pre-connected catheters (REF: 05436709, 05436717, 05436725). The Celsite® Dual Venous, Celsite® Brachial, and pediatric Celsite® Babyport® and Celsite® Babyport® S are included in the product line.
- Celsite® Access Ports (Dual Venous) are standard sized double access ports manufactured with titanium/epoxy/silicone materials that allow direct connection to a two-lumen silicone catheter. One catheter connection attaches both lumens of the catheter to the port (REF: 05430100, 05430101).
- Celsite® Access Ports (Brachial Venous) are small sized ports manufactured with titanium/epoxy/silicone materials and a polyurethane catheter. Celsite® brachial ports are for peripheral implantation in the arm of the patient (REF: 05433734).
- Celsite® *Discreet*™ Access Ports (Venous) are standard sized ports manufactured with a low-profile design with a left or right off-set exit cannula that permits a vertical incision and a 90 degree port/exit cannula. Celsite® *Discreet*™ ports are manufactured with titanium/epoxy/silicone materials. Catheters are available in either polyurethane or silicone materials (REF: 05430144, 05430145, 05430146).
- Celsite® *Concept*™ Access Ports (Venous) are standard sized ports manufactured with a low-profile, light weight design and have three silicone suture areas, which provide the physician with flexibility in suturing. Celsite® *Concept*™ ports are manufactured with titanium/polysulphone/silicone materials. Catheters are available in either polyurethane or silicone materials (REF: 05437025, 05437026, 05437028).
- Celsite® Babyport® and Babyport® S (Pediatric Venous) are small sized ports manufactured with titanium/epoxy/silicone materials. The Babyport® has a 4.5Fr polyurethane catheter; the Babyport® S has a 6Fr silicone catheter (REF: 05433742, 05433842).

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

### III – INDICATIONS

Celsite® Access Ports are indicated for patient therapies requiring repeated access to the vascular system.

- Celsite®/Celsite® **Concept™** Access Ports are indicated for intra-venous administration of drugs for chemotherapy, antibiotics and anti-viral drugs, and for parenteral nutrition, blood sampling or transfusions.
- Celsite® *Discovee™* / Celsite® Dual Venous / Celsite® Brachial / Celsite® Babyport® / Celsite® Babyport® S Ports are recommended for use whenever patient therapy requires repeated I.V. Access for injection, drug therapy and/or blood sampling.

### IV – CONTRAINDICATIONS

The following are contraindications to port placement:

- Known infection, bacteraemia or septicemia.
- Known allergy to any of the materials contained in the access port or catheter.
- If the medications to be used in the access port are incompatible with any of the materials contained in the access port or catheter.
- If the patient's anatomy does not allow the insertion of the catheter into the chosen access site.
- Blood clotting disorder.

### V – POTENTIAL COMPLICATIONS (Immediate and delayed)

Air embolism Brachial plexus injury Cardiac arrhythmia, tamponade, puncture Catheter disconnection or fragmentation Catheter occlusion Catheter rupture Embolism Explantation of the port secondary to inflammatory reaction Fibrin sheath formation Hematoma Hemothorax Hemorrhage	Infection/Sepsis Pneumothorax Port/Catheter migration Thromboembolism/Thrombophlebitis Thrombosis Fever, chills, or other flu-like symptoms Swelling, pain, pus, redness or burning sensation around the port incision Shortness of breath or dizziness Chest pain Pain or swelling of neck, face, or arm on side of body where the implantable port is inserted
--	---

Before implantation of the access port system these, and other, well documented potential complications should be considered.

### VI – GENERAL WARNINGS AND PRECAUTIONS

#### VI-1 HANDLING AND STORAGE

##### Precautions:

- This product should be implanted or accessed only by personnel experienced in the technical and clinical aspects of access ports.
- During storage avoid freezing or excessively high temperatures, storage at room temperature is recommended.
- Celsite® Access Ports are for single use only, do not re-sterilize the product, or any of the components, and destroy after use. The product is sterilized using Ethylene Oxide and is sterile and non-pyrogenic in unopened, undamaged, individual packages. Do not use if package is damaged.
- Check the expiration date before use and never implant a device past the expiration date.

#### VI-2 PATIENT

##### Precautions:

- Ensure that use of the access port is indicated taking into consideration the patients' anatomy, needs and proposed use of the port.
- Check the condition of the skin over the proposed implantation and puncture sites.
- For venous catheter placement, place the patient slightly inclined head down.
- Check hemostasis: Partial thrombin, Partial thrombin time, Bleeding time, Platelet count.

#### VI-3 CATHETER POSITION

**Precaution:** To ensure that the system works correctly, there should be no kinking of the catheter.

It is recommended that the catheter is tunnelled sub-cutaneously to the port.

#### VI-4 DRUG INCOMPATIBILITIES

**Precaution:** Do not use drugs or solvents which are incompatible with the raw material of the catheter.

Never use incompatible drugs together or in sequence: the port should be rinsed of all residues with 10-20 mL sterile sodium chloride (NaCl) 0.9% between each medication.

#### VI-5 NEEDLES

**Precautions:** Always use special bevelled needles which are adapted for use with implantable access port systems.

Normal hypodermic needles will damage the septum and may cause leakage of the system or blockage due to small silicone particles. Care must be taken to avoid needle sticks in accordance with local Infection Control standards (and avoid the risk of exposure to contaminated blood).

## VI-6 IMPLANTATION TECHNIQUES

### Precautions:

- To prevent shearing of the guidewire never withdraw the J guidewire through the Seldinger needle.
- Remove the guidewire and the dilator together. Do not remove the guidewire through the dilator as this may result in the guidewire unravelling.
- It is recommended that sterile NaCl 0.9% (with or without heparin according to local protocols) be used to flush the port and catheter during implantation (Fig. 3, page 14) to reduce the risk of thrombosis.
- During implantation ensure that the catheter is not damaged by unguarded forceps, suture needle or other sharp instruments.
- If using the Celsite® **Concept™** Access Ports (with silicone inserts) it is recommended that the suture be placed within the square indicated on the silicone to reduce the risk of port flip.

## VI-7 OCCLUSION

**Precaution:** Always use a syringe equal to or greater than 10 mL volume when rinsing the port to avoid generating excessive pressures in case of occlusion of the system.

**Warning:** In case of obstruction of the system never try to clear the blockage using a fluid under high pressure which carries the risk of catheter fracture and migration.

- According to local protocols, and under medical supervision, 70% alcohol may be used to aid the unblocking of silicone catheters when the blockage is due to lipid deposits. The use of alcohol with Polyurethane (PUR) catheters is not recommended.
- According to local protocols, and under medical supervision, Hydrochloric acid (HCl) 0.1 mol/L may be used to aid the unblocking of both silicone and PUR catheters when the blockage is due to mineral deposits.

## VI-8 INFECTION

**Warning:** In case of infection, if appropriate, treat with antibiotic drugs. If this fails, or is not appropriate, the catheter and port should be removed.

## VII CENTRAL VENOUS IMPLANTATION TECHNIQUES

### VII-1 Venous catheters

- Precaution:** The catheter tip should always be positioned in the Superior Vena Cava at the entrance to the right atrium (see Fig. 1, page 14). A radiograph should be taken to verify catheter position and exclude hemothorax or pneumothorax.
- Warning:** Particular attention should be paid when the catheter is to be implanted via the subclavian route. It is recommended that the catheter be inserted outside the costoclavicular space (see Fig. 2, page 14). Catheter ruptures have been observed with the subclavian route, the associated risk of extravasation of the infused drugs and embolisation of the distal extremity may have serious consequences. This is due to rupture of the catheter secondary to the catheter being pinched in the costoclavicular space (Pinch-Off syndrome). It is necessary to pay particular attention to catheters used for long periods or for ambulatory patients rather than short term catheters. Catheter rupture can occur with silicone and polyurethane catheters.

The following clinical signs may suggest catheter pinching:

- The patient needs to lift the arm to permit infusion.
- Intermittent malfunction of the catheter, such as difficulty with aspiration or infusion.
- Subclavicular pain or swelling during infusion.
- Palpitations or chest discomfort may indicate catheter fracture.

Remove any subclavian catheter which presents any of the signs of catheter Pinch-Off.

A radiograph with contrast injection may be useful to detect catheter abnormalities at the costoclavicular level.

### VII-2 Percutaneous Technique

- Insert the Seldinger needle into the chosen vein; verify the position by observing blood reflux (Fig. 4, page 14).
- Insert the J guidewire into the vein, when the correct position is reached, withdraw the needle (Fig. 5, page 14).
- Thread the assembled dilator and peelable sheath over the guidewire using a twisting motion to pass through the skin planes (Fig. 6, page 14).
- Remove the dilator and guidewire from the vein and insert the catheter through the peelable sheath (Fig. 7, page 14) to the desired position at the junction of the Superior Vena Cava and the right atrium (Fig. 1, page 14). Check catheter tip position using fluoroscopy.
- Prepare the port pocket at the chosen site, the port should lie approximately 1/2 - 1 cm below the skin surface away from the injection site. Cut the luer connector from the catheter and attach to the tunnelling rod (Fig. 8, page 14).
- Tunnel the catheter from the puncture point to the port pocket (Fig. 9, page 14); ensure there is no kinking of the catheter. Cut the excess catheter (at right-angles) prior to connecting to the port.
- Catheter port/connection:** Slide the connection ring over the catheter, firmly push the catheter onto the exit cannula ensuring the catheter covers the length of the exit cannula, slide the connection ring over the catheter and exit cannula. The connection ring should be in contact with the port (Fig. 10, page 14).
- Confirm catheter connection by gently pulling on the catheter.
- Insert the port into the port pocket, paying attention not to kink the catheter.
- Confirm catheter patency (by ensuring that both aspiration and injection are possible), taking care that the skin incision is not at the injection site.
- It is possible to fix the access port to the fascia with sutures.

### VII-3 Surgical Cut-down Technique

- The skin is incised and the tissue dissected (Fig. 11, page 14).
- The selected vessel is identified and tapes passed under it. The vessel is incised with a scalpel or with scissors (Fig. 12, page 14).
- The previously flushed catheter is inserted into the vessel, using either a small vessel dilator or vein lifter (Fig. 13, page 14).
- Follow from point e) - k) of the preceding paragraph: "Percutaneous Technique".

### VII-4 References requiring modifications to the above techniques

#### VII-4-1 Brachial Implantation Technique

Brachial ports and catheters can be placed using either the Percutaneous or Surgical Cut-down Technique, see above.

Select the appropriate vein on the anterior aspect of the elbow, normally the Basilic or Cephalic vein.

The catheter length should be measured before insertion of the catheter into the vein.

Ensure that the access port is not placed over the mid-line of the selected vein. The port and catheter should be positioned above the elbow on the internal surface of the arm.

#### VII-4-2 Braunule Technique

- Puncture the vein with a needle percutaneously or through a skin incision.
- Flush the catheter with sterile sodium chloride (NaCl) 0.9 %. Remove the stilette from the cannula, and insert the catheter through the cannula to the desired position at the junction of the Superior Vena Cava and the right atrium.
- Gently withdraw the cannula over the catheter.
- Follow from point e) - k) above: "Percutaneous Technique".

#### VII-4-3 Valved catheter

Do not place the catheter over the guidewire, as this may cause damage to the valve.

The catheter must be tunnelled from the puncture point to the port pocket; the tunnelling rod must not be placed through the valve.

#### VII-4-4 Pre-connected catheter and port

The catheter should be measured and trimmed to the correct length before insertion in the vein.

The catheter must be tunnelled from the port pocket to the puncture point.

## VIII – USE OF THE PORT AND CATHETER

### VENOUS PORTS

Always verify that the port and catheter are functional by aspirating a few mL of blood into a syringe and injecting a few mL of sterile sodium chloride (NaCl) 0.9% before attempting to start an infusion.

If aspiration of blood is not possible, attempt to inject a few mL of sterile sodium chloride (NaCl) 0.9% into the port.

**Caution:** If resistance to injection is noted, or if swelling occurs around the port or along the catheter, device malfunction should be suspected. In this case, device integrity should be verified using X-ray and contrast media studies.

## IX – MAINTENANCE OF THE PORT AND CATHETER

### General recommendations:

- If sterile heparinized sodium chloride (NaCl) 0.9% is used the system should be rinsed first with sterile sodium chloride (NaCl) 0.9% alone, as some drugs may react with heparin and result in blockage of the port/catheter due to the formation of precipitates.
- Always rinse the system between administration of different drugs. Special care should be taken with drugs which have a higher risk of precipitation, with anti-coagulation agents, after blood sampling or after transfusion to reduce the risk of catheter occlusion.
- Failure to maintain the system may result in occlusion of the catheter.

### IX-1 RINSING AND HEPARINIZATION

#### VENOUS PORTS

Rinsing of the access port is essential.

Rinse the port with 10-20 mL of sterile sodium chloride (NaCl) 0.9% after each use and every 4-6 weeks when no treatment is being given.

### IX-2 NEEDLE WITHDRAWAL

To avoid blood reflux into the catheter, always maintain positive pressure on the syringe during needle withdrawal (see Fig. 14, page 14).

## X – DURATION OF IMPLANTATION AND REMOVAL OF THE PORT

### X-1 DURATION OF IMPLANTATION OF THE SYSTEM

The port and catheter system should be removed at the end of treatment.

The duration of the system is mostly dependant on the catheter. Fracture is the most common complication of aging. Risk of fracture may increase over time depending on the material of the catheter, the pathway, and the entry point of the catheter into the vessel.

### X-2 REMOVAL OF THE SYSTEM

**Precaution** : When removing the system, care must be taken not to fracture the catheter. If the catheter is sutured into the vessel, the sutures should first be removed. Control the catheter while removing the system from the port pocket.

Be vigilant if there is excessive resistance to the removal of the catheter. Catheters may become encapsulated and attached to the vein wall. Should this occur and it is not possible to remove the catheter without risking catheter fracture, or if the catheter is fractured, the advice of an Interventional Radiologist, Surgeon or other Physician with endo-luminal experience should be sought.

### X-3 SPECIAL ATTENTION: BABYPORT® (pediatric), BABYPORT® S (pediatric) AND SMALL PORTS

Following implantation of a port in a child particular attention should be given to the position of the catheter tip. Over time, as the child grows, the position of the catheter tip will change and move higher in the Superior Vena Cava (SVC). As for any catheter, a high tip position in the SVC predisposes the catheter to fibrin sheath formation and other mechanical difficulties.

The position of the catheter tip should be checked radiographically a minimum of every 12 months or more frequently if the child is growing quickly. The catheter should be changed when the tip reaches the level of T4. Always ensure that the catheter is functioning correctly before beginning infusions.

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**B. Braun Interventional Systems Inc.**
en: Manufactured by : **B. BRAUN MEDICAL SAS**