

ACCEL™ Valved Safety Centesis Catheter

with Introcan Safety™ Technology



The Valved Safety Centesis Catheter will be used to drain fluid during thoracentesis and paracentesis procedures. There are three drainage holes located around the catheter to maximize the drainage of fluid. The catheter surface and tapered end is designed for smooth transition through the tissue, while the inner surface is designed for easy retraction of the insertion needle. The insertion needle comes with a vent plug that allows air to exit and does not allow fluid to exit, providing safety for the clinician. The echogenic needle tip provides guidance for proper catheter placement.

Indications for Use

- The Valved Safety Centesis Catheter is indicated for the therapeutic or diagnostic aspiration of fluids from body cavities.

The Valved Safety Centesis Catheter is designed to:

- Minimize inadvertent needlestick injuries with a passive safety needle-shielding feature.
- Control the flow of fluid, which aids in the prevention of fluid exposure.

STERILE EO in unopened, undamaged package.
NONPYROGENIC Rx only



For single use only. Not made with natural rubber latex or DEHP.
Do not resterilize.

Pictorial is for reference only.

A4807754-1 ACCEL is a trademark of B. Braun Interventional Systems Inc.
REV. 1/16

B | BRAUN

B. Braun Medical Inc.

Bethlehem, PA 18018-3524 USA
1-800-227-2862

NOTES:

- Only persons knowledgeable of the risks involved and qualified in the procedures performed should use this device.
- The device is luer compatible.
- Report needle stick injuries immediately, and follow established protocols.

CAUTIONS:

- IN THE UNLIKELY EVENT THAT THE SAFETY SHIELD DOES NOT ENGAGE, KEEP THE NEEDLE POINT AWAY FROM THE BODY AND FINGERS AT ALL TIMES, AND IMMEDIATELY DISPOSE INTO AN APPROVED SHARPS CONTAINER.
- NEVER REINSERT THE NEEDLE INSIDE THE CATHETER ONCE THE NEEDLE HAS BEEN PARTIALLY OR COMPLETELY WITHDRAWN, AS IT MAY PIERCE AND/OR SEVER THE CATHETER.
- DO NOT ATTEMPT TO OVERRIDE OR DEFEAT THE SAFETY SHIELD.

Instructions For Use: Use Aseptic Technique.

1. Raise skin wheal.
2. Remove needle guard.
3. Perform puncture using the Valved Safety Centesis Catheter.
4. Introduce the Valved Safety Centesis Catheter into the cavity.
5. If applicable, confirm placement of needle tip using ultrasound guidance or fluoroscopy.
6. Verify flashback of fluid into the needle hub.
7. If applicable, reconfirm placement of needle using ultrasound guidance or fluoroscopy.
8. Remove insertion needle from catheter.
9. Connect to drainage collection device.
10. After drainage procedure is complete, remove catheter and discard.