B. Braun Inflation Device

GB B. Braun Inflation Device

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GB

INDICATION

This Inflation Device is recommended for use during the performance of balloon angioplasty procedures to inflate and deflate the balloon and to measure the pressure within the balloon during the procedure.

DESCRIPTION

The disposable B. Braun Inflation Device consists of a 25 cc syringe barrel with a pressure gauge, a threaded plunger assembly with Winged Lock Mechanism used to generate and control balloon inflation pressures, a flexible high pressure extension tube, and three-way high pressure stopcock. The gauge is calibrated from 0 to 30 atm (0 to 441 psi) of positive pressure. The accuracy of the gauge is within 1 atm over the range. Accuracy of syringe graduations: ±5%

WARNING

Isovue® 370 or Renovist™ II in combination with saline may cause copper oxide to enter solution. Care should be used with copper-sensitive patients. Refer to the directions accompanying the balloon dilation catheter for instructions for use, maximum pressure, precautions, and warnings for that device.

INSTRUCTIONS FOR USE: Use Aseptic Technique.

Inspect prior to use to insure no damage has occurred during shipping and handling.

NOTE: Product is packaged in a sterile peel pouch to facilitate sterile transfer of the device.

Inspect for damage prior to use. Follow steps below to check inflation device.

1. Turn winged lock mechanism counter-clockwise so wings are parallel to the gauge.
2. Aspirate air by pulling plunger back to 25cc graduation.
3. Attach stopcock to extension tube of inflation device and turn stopcock handle with the "OFF" position facing extension tubing port.
4. Lock syringe by turning the winged lock mechanism clockwise 90° so wings are perpendicular to the gauge.
5. Turn the plunger handle clockwise and advance the syringe plunger forward to build pressure.
6. Pressurize inflation device to approximately 44.1 psi (3 atm), let gauge stabilize 10 seconds and observe gauge for any loss of pressure.
7. Do not use if damage or loss of pressure is observed.

PREPARATION

1. Prepare solution of contrast medium and normal saline as recommended by catheter manufacturer.
2. Insert the end of the extension tube into the contrast solution.
3. Aspirate up to 20 cc of solution into the Inflation Syringe.
4. Purge all air from the inflation syringe barrel by pointing it upwards so that the air accumulates in the distal end of the barrel. Slowly advance the syringe plunger until all air has been expelled from the syringe barrel and extension tube. CAUTION: Inspect the syringe tubing and stopcock (if used) to insure that there is no air in the system. Recommended optimal volume of contrast solution is 8 to 10 cc.

ATTACHING THE INFLATION DEVICE TO THE BALLOON

1. Prepare and test the balloon catheter according to the manufacturer's directions for use.
2. Create a fluid-fluid connection between the balloon and the Inflation Syringe extension tube and connect the luer connectors securely.
3. Unlock by turning the Winged Lock Mechanism 90° counter-clockwise and pull back on the plunger handle to apply a vacuum to the balloon.

INFLATION INSTRUCTIONS

1. To engage, turn the winged lock mechanism clockwise 90°. Ensure the latch is fully engaged.
2. Slowly rotate the plunger clockwise to increase the pressure or counterclockwise to decrease the pressure. WARNING: A slight rotation of the handle will change the pressure significantly. Monitor pressures continually to avoid overpressurization.

DEFATION INSTRUCTIONS

1. Rotate the winged lock mechanism counterclockwise 90° to a horizontal position. Pressure will be quickly released as the plunger is fully withdrawn and all fluid is removed from the balloon. NOTE: Negative pressure can be maintained by engaging the winged lock mechanism when the plunger is fully withdrawn. CAUTION: Observe the analog gauge to confirm that the system maintains a vacuum. Loss of negative pressure indicates the presence of a leak in the system.

PRECAUTIONS

1. Follow all instructions recommended by the balloon dilation catheter manufacturer.
2. Use the sterile contrast media that the balloon dilation catheter manufacturer recommends.
3. The balloon dilation catheter manufacturer's recommended maximum balloon inflation pressure should not be exceeded.
4. DO NOT exceed a pressure of 441 psi (30 atm) when using this device. WARNING: Refer to the manufacturer's directions accompanying the balloon dilation catheter for specific information on use, maximum inflation pressure, precautions, and warnings for that device. Inflation pressures should be closely monitored when inflating balloon. The inflation syringe is a high volume, low compliance system capable of generating high pressures with relative ease.

US PATENT NO.: 5,860,955

ISOVUE is a registered trademark of Bracco Diagnostics, Inc.
RENOVIST is a trademark of Bracco Diagnostics, Inc.

CAUTION: Read the Instructions for Use of all equipment to be used in the procedure prior to use.

WARNING

Re-use of single-use devices creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/or limited functionality of the device may lead to injury, illness or death of the patient.

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