---

## Multi-Snare® Multi-Snare® Micro Multi-Snare® VSD Loop-Set

---

en	Instructions for Use	5
de	Gebrauchsanweisung	8
fr	Mode d'emploi	11
es	Instrucciones de uso	14
it	Istruzioni per l'uso	17
nl	Gebruiksaanwijzing	20
pt	Instruções de utilização	23
sv	Bruksanvisning	26
el	Οδηγίες χρήσεως	29
tr	Kullanma kılavuzu	32
cs	Návod k použití	35
hu	Használati útmutató	38
ru	Инструкции	41
pl	Instrukcja użytkowania	45
ar	إرشادات الاستخدام	48

---

Explanation of the symbols on label and packaging // Erklärung der Symbole auf Etikett und Verpackung // Explication des symboles sur l'étiquette et l'emballage // Explicación de los símbolos en la etiqueta y el envoltorio // Illustrazione dei simboli riportati su etichetta e confezione // Verklaring van de symbolen op het etiket en de verpakking // Explicação dos símbolos na etiqueta e embalagem // Förklaring av symboler på etiketten och förpackningen // Εξήγηση των συμβόλων πάνω στην ετικέτα και τη συσκευασία // Etiket ve ambalaj üzerindeki sembollerin açıklamaları // Vysvětlivky k symbolům na etiketě a balení // A címkép és a csomagoláson lévő szimbólumok magyarázata // Объяснение символов на этикетке и упаковке // Objaśnienie symboli na etykietcie i opakowaniu. // شرح الرموز الموجودة على الملصق والعبوة

**REF**

Reference number // Bestellnummer // N° de commande // Número de referencia para pedidos // Numero d'ordine // Artikelnummer // Número de referência // Bokningsnummer // Αριθμός παραγγελίας // Sipariş numarası // Objednací číslo // Rendelési szám // Номер заказа // Numer zamówienia // الرقم الاستدلالي

**LOT**

Lot number // Chargennummer // N° de lot // Número de lote // Numero di lotto // Chargennummer // Número de lote // Partinummer // Αριθμός παρτίδας // Parti numarası // Číslo šarže // Gyártási tétel száma // Номер партии // Numer partii // رقم التشغيل

Read Instructions for Use carefully // Gebrauchsanweisung beachten // Respecter les instructions du mode d'emploi // Lea atentamente las instrucciones de uso // Osservare le istruzioni per l'uso // Volg de instructies in de gebruiksaanwijzing // Leia atentamente as instruções de utilização // Läs bruksanvisningarna noggrant // Διαβάστε προσεκτικά τις οδηγίες χρήσεως // Kullanna kilavuzunu dikkate alınız // Respektujte Návod k použití // Vegye figyelembe a használati útmutatót // Соблюдать инструкцию по применению // Przestrzegać instrukcji użytkowania // اقرأ إرشادات الاستخدام بعناية

Protect from direct sunlight // Vor direkter Sonneneinstrahlung schützen // Protéger de l'exposition directe au soleil // Proteger de la luz solar directa // Proteggere dall'irradiazione solare diretta // Product niet blootstellen aan direct zonlicht // Proteger da luz direta do sol // Skydda från direkt solljus // Να μην εκτεθούν στο άμεσο ηλιακό φως // Doğrudan güneş ışığından koruyunuz // Chraňte před přímým slunečním zářením // Όντα η közvetlen napfénytől // Не допускать попадания прямых солнечных лучей // Chronić przed bezpośrednim oddziaływaniem promieni słonecznych // حافظ على حماية المنتج من أشعة الشمس

Store in a dry place // Trocken aufbewahren // A conserver au sec // Guardar en lugar seco // Conservare in un luogo asciutto // Droog bewaren // Armazenaar num local seco // Förvara torr // Να διατηρείται σε ξηρό μέρος // Kuru ortamda saklayınız // Uchovávejte v suchu // Száraz helyen tartandó // Хранить в сухом месте // Przechowywać w suchych pomieszczeniach // التخزين في مكان جاف

Expiry Date // Verfalldatum // Date de péremption // Fecha de caducidad // Data di scadenza // Vervaldatum // Prazo de validade // Sista förbrukningsdag // Ημερομηνία λήξεως // Son kullanna tarihi // Použitelné do // Lejárati idő // Срок хранения // Data ważności // تاريخ انتهاء الصلاحية

For single use only // Nur zum einmaligen Gebrauch bestimmt // à usage unique seulement // Para un solo uso // Solo per impiego monouso // Alleen bestemd voor eenmalig gebruik // Utilização única // Använd endast en gång // Για μία χρήση μόνο // Bir defahk kullanim içindir // Určeno k jednorázovému použití // Kizárólag egyszeri használatra alkalmas // Только для одноразового использования // Przeznaczone do jednorazowego użytku // للاستخدام مرة واحدة فقط

Do not resterilize // Nicht erneut sterilisieren // Ne pas restériliser // No volver a esterilizar // Non risterilizzare // Niet opnieuw steriliseren // Não reesterilizar // Återsterilisera inte // Μην αποστειρώνετε εκ νέου // Yeniden sterilize edilemez // Nesterilizujte znovu // Nem újrasztirilizálható // Не подвергать рестерилизации // Nie poddawać ponownej sterylizacji // لا تقم بإعادة التعقيم

**STERILE EO**

Sterilized by Ethylene Oxide // Sterilisiert mit Ethylenoxid // Stérilisé à l'oxyde d'éthylène // Esterilizado con óxido de etileno // Sterilizzato con ossido di etilene // Gesteriliseerd met ethyleenoxide // Esterilizado por óxido de etileno // Steriliserat med etylenoxid // αποστειρωμένο με αιθυλενοξειδίο // Etilen oksit ile sterilize edilmiştir // Sterilizovano etylenoxidem // Etílenoxid-dal sterilizált // Простерилизовано оксидом этилена // Sterylizowane tlenkiem etylenu // معقم بإكسيد الإيثيلين

Do not use if packaging is damaged // Bei beschädigter Verpackung nicht verwenden // Ne pas utiliser si l'emballage est endommagé // No usarlo si el embalaje está dañado // Non utilizzare se la confezione risulta danneggiata // Niet gebruiken indien de verpakking beschadigd is // Nunca utilize se a embalagem estiver danificada // Använd inte om förpackningen är skadad // Μην το χρησιμοποιείτε εάν η συσκευασία είναι χαλασμένη // Ambalajın hasarlı olması halinde kullanmayınız // Nepoužívejte v případě poškození obalu // Sérült csomagolás esetén ne használja a terméket // Не использовать, если упаковка повреждена // Nie stosować w przypadku uszkodzenia opakowania // لا تستخدم المنتج إذا كانت العبوة تالفة

Manufacturer // Hersteller // Fabricant // Fabricante // Produttore // Fabrikant // Fabricante // Tillverkare // Κατασκευαστής // Üretici // Výrobce // Gyártó // Изготовитель // Producent // جهة التصنيع

Latex free // Latex-frei // Sans latex // Sin látex // Non contiene lattice // Latex-vrij // Isento de látex // latexfri // Χωρίς λάτεξ // Lateks içermez // Neobsahuje latex // Latexmentes // Не содержит латекса // Nie zawiera lateksu // خال من اللاتكس

Diethylhexylphthalate (DEHP) free // Diethylhexylphthalate-(DEHP)-frei // Sans diéthylhexylphthalate (DEHP) // Sin dietilhexilftalato (DEHP) // Non contiene dietilesilftalato (DEHP) // Zonder diethylhexylftalato (DEHP) // Isento de dietilhexilftalato (DEHP) // fri från diethylhexylftalát (DEHP) // Χωρίς δι(2-εθυλεξυλ) φθαλικό εστέρα (DEHP) // Dietil heksilftalát (DEHP) içermez // Neobsahuje diethylhexylftalát (DEHP) // Dietilén-ftalátot (DEHP) tartalmaz // Не содержит диэтилгексилфталат (ДЭГФ) // Nie zawiera DEHP // (DEHP) خال من ثنائي إيثيل هيكسيل الفثالات

Rx only

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician. // Hinweis: Laut Bundesgesetzgebung der USA ist die Verwendung dieses Produkts nur durch einen Arzt gestattet // Remarque : Selon la législation fédérale américaine, l'utilisation de ce produit est autorisée uniquement par un médecin // Nota: La normativa federal vigente (EE.UU.) limita el uso de este dispositivo a médicos // Nota: La legge federale (USA) limita l'uso di questo dispositivo a medic // ) // Opmerking: Volgens de federale wet in de VS mag dit product uitsluitend door een arts worden gebruikt. // Aviso: Conforme a legislação federal dos EUA, a utilização deste produto está apenas autorizada a médicos // Varning: Federal lag (USA) tillåter endast användning av denna produkt för läkare // Προσοχή: Σύμφωνα με την ομοσπονδιακή νομοθεσία των ΗΠΑ η χρήση αυτού του προϊόντος είναι επιτρεπτή μόνο από έναν ιατρό // Upozornění: Podle národní legislativy USA je dovoleno používat tento produkt pouze lékaři // Tájékoztató: Az Amerikai Egyesült Államok szövetségi törvénye alapján a termék alkalmazása csak orvos által engedélyezett // Указание: Согласно федеральному законодательству США, данным продуктом может пользоваться только врач // Informacja: Zgodnie z prawem federalnym Stanów Zjednoczonych do stosowania niniejszego produktu uprawnieni są wyłącznie lekarze. // تحذير: يحظر القانون الفيدرالي (بالولايات المتحدة الأمريكية) بيع هذا الجهاز إلا من خلال الطبيب أو بناء على طلبه.

## Instructions for Use

### Multi-Snare<sup>®</sup>, Multi-Snare<sup>®</sup> Micro, Multi-Snare<sup>®</sup> VSD Loop-Set

#### Product description

The Multi-Snare<sup>®</sup> and the Multi-Snare<sup>®</sup> Micro are both considered as medical devices. They merely differ because of their dimensions. They are made out of a super elastic Nitinol wire, wrapped with a thin, radiopaque Platinum wire around the distal pre-formed loop. A PTFE tube is shrunk on the proximal shaft of this loop to ensure good torqueability. The super elastic characteristics and the shape memory effect of the Nitinol material allow the snare to be inserted through the pfm medical snare catheter or through another catheter (check compatibility prior use) to the region of interest without any risk of deformation. Depending on the position of the snare catheter's tip the diameter of the loop is variable in a predetermined range (see label).

Both Multi-Snare<sup>®</sup>, and Multi-Snare<sup>®</sup> Micro are available as a set including the snare catheter. For better radiopacity a radiopaque marker is located distally on the snare catheter shaft.

An insertion aid is provided to ease the guidance of the snare into the catheter. The provided torquer eases the handling of the product.

The Multi-Snare<sup>®</sup> Micro is available with diameters of 2 and 4 mm and in the length of 175 cm, whereas the standard Multi-Snare<sup>®</sup> is offered from 5 to 30 mm in diameter and in the length of 125 cm.

The Multi-Snare<sup>®</sup> VSD Loop-Set contains a 15 mm snare as described above, but with an extended shaft. Also, the set features a 260 cm hydrophilic-coated crossover wire (0.035 inch with 60° angled tip), and 2 metal clips.

**The Multi-Snare<sup>®</sup> VSD Loop-Set is not for use in the USA.**

Please refer to the product label for exact dimensions, designs and compatibility with other devices.

#### General Information and Warnings

The product described here may be used only by specialist physicians who are familiar with the techniques used in diagnostic and interventional cardiology.

These instructions for use and the information on the packaging should be read carefully before each use.

#### WARNINGS

-  **Federal Law (USA) restricts this device to use by or on the order of a physician.**
-  **This product contains chemicals known to the State of California to cause cancer and birth defects or other reproductive harm.**
-  **Do not use the product if the packaging has been opened, or is damaged, if you are not sure that it is sterile, or if the expiry date has passed.**
-  **Each product is packed separately, and is delivered in an EO-sterilised and pyrogen-free condition. It is intended for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization of medical devices intended for single use only may result in degraded performance or a loss of functionality. Reuse of single use only medical devices may result in exposure to viral, bacterial, fungal or prionic pathogens.**
-  **The product must be stored in a dry place. Do not expose the packaged products to direct sunlight.**
-  **Do not clip the Multi-Snare<sup>®</sup> VSD Loop-Set's clips directly onto the wire. Always protect the wire with compresses against direct contact with the clips in order to avoid damaging the hydrophilic coating.**

#### Indications

- In general for the removal of foreign bodies out of the vascular system or a hollow organ, and for the manipulation of catheters, guide wires etc. in the vascular system or in a hollow organ.
- Because of its dimensions (small diameter and long shaft) the Multi-Snare<sup>®</sup> Micro is especially suited for the removal of foreign bodies out of distal smaller vessels.
- The Multi-Snare<sup>®</sup> VSD Loop-Set is especially intended for supporting the contra-lateral positioning of a catheter using the crossover technique by aligning the guide wire. The shaft of the snare is extended especially for this indication in order to be able to push it forward over a lying flow-directed balloon-tipped catheter. The clips are suitable for fixing the wire in the crossover position.

Other examples of manipulation using a snare:

- Assistance at veni-puncture procedure in order to provide central venous access.
- Fibrin sheath stripping of an indwelling venous catheter.
- Positioning resp. correction of position of an indwelling venous catheter.

### Indications (USA and Canada Only)

- Retrieval and manipulation of foreign objects from the vascular system and hollow viscera.
- Assistance in creating loops where cross-over technique is applied.
- Repositioning of indwelling venous catheters.
- Assistance in performing venipuncture to obtain access to central vein.

### Contraindications

The snare may not be used

- For removal of foreign bodies that have become entrapped by tissue overgrowth.
- For removal of implanted pacing leads.
- For fibrin sheath stripping in the presence of septal defects of persistent Foramen Oval (Reason: Risk of Cross embolism).

### Possible Complications

- Air embolism
- Death
- Haematoma at the puncture site
- Infection
- Perforation of the vessel wall
- Pulmonary embolism because of “lost” fibrin sheath material, especially when the fibrin sheath is long, and big

### Directions for use

#### Selection of the snare, and preparation

The maximum snare diameter (see label) should not exceed the targeted vessel diameter.

Remove the snare from the protective transport dispenser. Take care that the proximal end of the snare does not cross the mark on the connected insertion aid.

Carefully inspect snare (and snare catheter) for any damages.

Note: It is recommended to use the pfm medical Multi-Snare® or Multi-Snare® Micro together with a pfm medical snare catheter.

Note: Prior to use thoroughly flush the snare catheter, and the Multi-Snare® / Multi-Snare® Micro with sterile heparinized saline solution. An additional systemic application of heparin is recommended.

#### Introduction

Disconnect the insertion aid from the transport dispenser, and connect it with the snare catheter or another suitable catheter.

**Note:** We advise you that the insertion aid was specially optimized for the compatibility with the pfm medical snare catheter.

#### Insertion without guide wire

The snare catheter is inserted into the human body via an introducer sheath or a big vascular cannula with haemostatic valve.

Insert the snare using the insertion aid provided into the snare catheter which is already positioned in the vascular system, and place it proximally to the foreign body. For an eased manipulation and torqueability of the snare the provided torquer should be mounted at the proximal end of the Multi-Snare®.

In general it is recommended to insert first the snare catheter to the region of interest. It is also possible to introduce the snare catheter with the preloaded snare as a unit.

#### Insertion over a guide wire

**Note:** Check whether the guide wire or guiding catheter have the necessary diameter, and length for catheter removal resp. to accommodate the snare. If necessary, replace with a suitable type.

Position the indwelling guide wire as close as possible to the foreign body. Place the snare loop over the proximal end of the guide wire, and pull tight. Advance the snare into the guiding catheter or vascular cannula until the distal end of the snare catheter is positioned just proximal to the foreign body.

#### Manipulation/Extraction of foreign objects

- ➊ Push the snare shaft carefully forward to open the loop. The open loop is then slowly advanced forward and placed around the proximal end of the foreign body. The diameter of the snare can be adjusted within the determined limits by moving the snare catheter forward or backward.
- ➋ By advancing the snare catheter, the snare loop is closed around the foreign body, and the foreign body is trapped.

**WARNINGS****⚠ Pulling back the snare can change the position of the snare loop.**

- ① To manipulate a foreign body, maintain tension on the snare catheter to retain the hold on the foreign body. Move the snare, and the snare catheter as a unit, in order to manipulate the foreign body to the desired position.
- ② To remove the foreign body, maintain tension on the snare catheter, and pull the snare and the snare catheter as a unit to the tip of resp. into the guiding catheter or introducer sheath. The foreign body is then withdrawn through or together with the guiding catheter or introducer sheath. The withdrawal of larger foreign bodies may require the exchange with larger sheaths, or guiding catheters, or a comminution of the foreign body in the peripheral vascular system.

**Snare Assisted Venous Canalization**

- ① Introduce the snare at a patent venous access site, and position it at the new, desired venous access site in the vascular system.
- ② Open the snare loop. The opened loop functions as a target to guide the puncture at the desired venous access site.
- ③ Insert the guide wire through the cannula, and the opened snare loop.
- ④ Pull the snare tight around the cannula by advancing the snare catheter.
- ⑤ Pull the guide wire into the inferior vena cava.

**Snare Assisted Fibrin Sheath Stripping**

- ① Advance the snare to the desired place from a common femoral vein access through the Vena Femoralis.
- ② Insert a suitable guide wire through the indwelling catheter, and forward it, until it exits the distal end.
- ③ Encircle the guide wire with the snare loop.
- ④ Advance the snare over the distal end of the catheter to a position proximal to the fibrin sheath.
- ⑤ Tighten the snare around the catheter. Hold the diameter of the snare loop by continuing to apply light traction while gently pulling the snare down (towards the distal end of the catheter) over the hub.
- ⑥ Repeat steps 4, and 5 until the catheter is free of fibrin sheath.

**WARNINGS**

- ⚠ **When using the snare for fibrin sheath stripping the user must ensure, that no high pulling force is applied to the catheter. Too high pulling force may lead to the damage and failure of the product, e.g. stretching, breakage, etc.**

**Disposal after Use**

After use, medical products and accessories pose a potential biological hazard. For this reason, the products and their accessories should be handled and disposed of in accordance with recognised medical procedure, and in compliance with the relevant legal regulations and local ordinances.

**Warranty**

Improper handling or use may damage the pfm medical Multi-Snare®. pfm medical disclaims all warranties, expressed and implied, including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. pfm medical shall not be liable to any person or entity for any medical expenses or any direct, incidental or consequential damages caused by any defect, failure or malfunction of the system, whether a claim for such damage is based upon warranty, contract, tort or otherwise.

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### Multi-Snare<sup>®</sup>, Multi-Snare<sup>®</sup> Micro, Multi-Snare<sup>®</sup> VSD Loop-Set

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