

## B. Braun Interventional Systems

### Celsite® *Discreet*™ Access Ports

Implantable Access Port Systems (venous)



Celsite® *Discreet*™ Access Ports for  
mid-to long-term venous access



- The unique anatomic delta-shape design of the Celsite® *Discreet*™ reduces scarring and maximizes patient comfort
- The 90° off-set port exit cannula ensures port stability

# Celsite® *Discreet*™

- Rapid, secure and easy implantation
- Lightweight anatomic design for easy insertion and patient comfort
- Excellent port stability due to the delta shape port design
- Radiopaque
- 90° port/exit cannula configuration, virtually eliminates the risk of port flippage and the need to suture
- Latex and DEHP-free
- Titanium chamber for enhanced security
- Catheters marked every cm from 5 cm

Size	Exit Cannula	Description	OD	OD (mm)	ID (mm)	Length (mm)	Flow Rates (mL/min)*		Model Number	Reference	Accessories Included
							19G	22G			
 <p>Standard</p>	Right Exit	Large Silicone	8.5F	2.8	1.1	800	26	11	STR201L	05430145	X
	Left Exit	Large Silicone	8.5F	2.8	1.1	800	26	11	STL201L	05430144	X
 <p>Small</p>	Left Exit	Small Silicone	6.5F	2.2	1.1	800	22	11	STL205F	05430146	X

\*Flow rates established by gravity infusion of NaCl 0.9%, height 1m, catheter length 40 cm.

## Ordering Information:

For more information or to place an order, contact your B. Braun Interventional Systems representative or call 1-877-VENA CAV (836-2228)

### Indications For Use:

Celcite® Access Ports are indicated for patient therapies requiring repeated access to the vascular system. Celcite® *Discreet*™ Access Ports are recommended for use whenever patient therapy requires repeated I.V. Access for injection, drug therapy and/or blood sampling.

Refer to the Celcite® *Discreet*™ Instructions for Use for relevant warnings, precautions, complications and contraindications.



- 18G thin wall needle (Seldinger)
- J guidewire
- Peelable introducer
- Vein lifter
- Tunneling rod
- 10 mL Omnifix® syringe
- 22G Surecan® needle (x2)
- 20G x 20 mm winged Surecan® needle
- Catheter
- Connection rings (x2)



DEHP free



Components and packaging do not contain natural rubber latex

Rx only

CV1082 1/14

©2014 B. Braun Interventional Systems Inc.



Manufactured by:  
B. Braun Medical (France)  
204 Avenue de Marechal Juin  
F-92100 Boulogne Billancourt  
FRANCE



Manufactured for  
B. Braun Interventional Systems Inc.  
824 Twelfth Avenue  
Bethlehem, PA 18018  
Tel: 1-877-VENA CAV (836-2228)  
Fax: 610-849-1334  
www.bisusa.org