

# DISCOVER NEW HORIZONS IN FLUID DRAINAGE

- The ASEPT® 600ml & 1,000ml EVACUATED drainage bottles are manufactured with NON-GLASS, shatterproof, medical grade plastic for safe and secure containment of bio-hazardous fluids.
- Each ASEPT® bottle is crafted for accuracy.
- Available with valved, Y-site drain line or without drain line for direct, needle-free closed system connection.







# ASEPT® 600ml & 1,000ml EVACUATED Drainage Bottle with Y-site Sampling Port and Drain Line

# ASEPT® 600ml & 1,000ml EVACUATED Drainage Bottle without Drain Line Available with or without tyvek pouch

MADE OF  
MEDICAL GRADE  
PLASTIC



Drainage Bottle without tyvek pouch (sterile fluid pathway only)

Drainage Bottle with tyvek pouch (sterile)

## Ordering Information:

ITEM	Ref Number	Case Qty.
ASEPT® 600ml Drainage Bottle with Y-site Sampling Port and Drain Line	M7052	10
ASEPT® 1,000ml Drainage Bottle with Y-site Sampling Port and Drain Line	M7050	10

**ASEPT® Drainage Bottle with or without Drain Line Indications for Use:**  
The ASEPT® 600ml or the ASEPT® 1,000ml Evacuated Drainage Bottle is a non-powered, single patient, portable suction device that consists of a manually operated plastic disposable evacuation system intended to provide a vacuum for suction drainage of surgical wounds. This set is used for large wound-areas. (i.e., after surgical operations)

**ASEPT® Bare Drainage Bottle Indications for Use:**  
The ASEPT® Drainage Bottle is a non-powered, single patient, portable suction device used to aspirate, remove, or sample body fluids.

Refer to the ASEPT® 600ml & 1,000ml Drainage Bottle Instructions for Use for relevant warnings, precautions, complications and contraindications. This device has been designed for single use only.

**Rx only**  
CV2005-1/14

ASEPT® is a registered trademark of pfm medical, inc.  
©2014 B. Braun Interventional Systems Inc.

## Ordering Information:

ITEM	Ref Number	Case Qty.
ASEPT® 600ml Evacuated Drainage Bottle without Drain Line	622272	10
ASEPT® 1,000ml Evacuated Drainage Bottle without Drain Line	622273	10
ASEPT® 600ml Bare Evacuated Drainage Bottle without tyvek pouch (sterile fluid path only)	622272B	10
ASEPT® 1,000ml Bare Evacuated Drainage Bottle without tyvek pouch (sterile fluid path only)	622273B	10

For more information or to place an order, contact your B. Braun Interventional Systems Inc. representative or call 1-877-VENA-CAV (836-2228)



Manufacturer:  
pfm medical, inc.  
1815 Aston Ave.,  
Suite 106  
Carlsbad, CA 92008

Manufactured for:  
B. Braun Interventional Systems Inc.  
824 Twelfth Avenue  
Bethlehem, PA 18018 USA  
Tel: 1-877-VENA CAV (836-2228) (USA)  
Fax: 1-610-849-1334  
www.bisusa.org