

Elite HV[®]


Hemostasis Valve Introducer System

Confidence and Selection From The Start



- Easy catheter movement through the hemostasis valve during multiple catheter insertions
- Kink-resistant material
- Smooth dilator-to-sheath transition for atraumatic access
- Available in 3-9F sizes in 7, 12 and 24cm lengths
- Radiopaque marker band available in 5-7F sizes
- DEHP and latex-free

Elite HV[®] Hemostasis Valve Introducer System



Interference Fit ensures a snug, smooth dilator-to-sheath transition, minimizing the potential for tissue trauma during vascular access

Snap-Fit Dilator Hub locks dilator and sheath together during advancement

Silicone gasket enables excellent device movement

Kink-resistant material—flexible design maintains sheath integrity

Radiopaque Marker Band to facilitate sheath insertion

Ordering Information:

REF Number	French Size	Sheath Length (cm)	Guidewire Shaft Material	Guidewire Tip Material	Guidewire Diameter & Length	Case Qty.
7cm length without Marker Band DUAL TIP						
612930	3	7	NT	PT/SS*	.018 x 40cm	10
612800	3	7	NT	W/SS*	.018 x 40cm	10
7cm length without Marker Band						
612801	3	7	SS	SS	.021 x 45cm	10
612802	4	7	SS	SS	.021 x 45cm	10
612803	5	7	SS	SS	.021 x 45cm	10
612804	6	7	SS	SS	.021 x 45cm	10
612805	7	7	SS	SS	.021 x 45cm	10
12cm length without Marker Band						
612806	4	12	SS	SS	.035 x 50cm	10
612807	5	12	SS	SS	.038 x 50cm	10
612808	6	12	SS	SS	.038 x 50cm	10
612809	7	12	SS	SS	.038 x 50cm	10
612810	8	12	SS	SS	.038 x 50cm	10
612811	9	12	SS	SS	.038 x 50cm	10
24cm length without Marker Band						Recommended Guidewire Dia.
612815	6	24	—	—	.038	5
612816	7	24	—	—	.038	5
612817	8	24	—	—	.038	5
612818	9	24	—	—	.038	5

* This kit contains a dual tipped wire

The Elite HV[®] hemostasis valve introducer system is indicated for use in percutaneous procedures to introduce catheters and other intravascular devices into the vasculature. Refer to the Elite HV[®] instructions for use for relevant warnings, precautions, complications and contraindications.

REF Number	French Size	Sheath Length (cm)	Guidewire Shaft Material	Guidewire Tip Material	Guidewire Diameter & Length	Case Qty.
7cm length with Marker Band						
612903	5	7	SS	SS	.021 x 45cm	10
612904	6	7	SS	SS	.021 x 45cm	10
612905	7	7	SS	SS	.021 x 45cm	10
12cm length with Marker Band						
612907	5	12	SS	SS	.038 x 50cm	10
612908	6	12	SS	SS	.038 x 50cm	10
612909	7	12	SS	SS	.038 x 50cm	10
24cm length with Marker Band						Recommended Guidewire Dia.
612915	6	24	—	—	.038	5
612916	7	24	—	—	.038	5

Guidewire Material Legend:

NT = Nitinol
W = Tungsten
SS = Stainless Steel
PT = Platinum

For more information or to place an order, contact your B. Braun Interventional Systems Inc. representative or call 1-877-VENA-CAV (836-2228)

Rx only
CV330- 7/12

©2012 B. Braun Interventional Systems Inc.



Components and packaging do not contain natural rubber latex



Manufacturer:
GALT MEDICAL CORP.
2220 Merritt Drive
Garland, TX 75041

B. BRAUN

Interventional Systems

Manufactured for:
B. Braun Interventional Systems Inc.
824 Twelfth Avenue
Bethlehem, PA 18018 USA
Tel: 1-877-VENA CAV (836-2228) (USA)
Fax: 1-610-849-1334
www.bisusa.org