

Intradyn™

Hemostasis Introducer

Confidence and selection
from the start.



Design Features:


- Easy catheter movement through the hemostasis valve during multiple catheter insertions.
- Leak-free hemostasis valve, with or without device in place.
- Kink resistant sheath.
- Smooth dilator to sheath transition for atraumatic access.
- Complete range of sizes and lengths – 10 through 14 French in 10 or 23 cm lengths.



Intradyn™ Hemostasis Introducer

Snap-Fit Dilator Hub locks dilator and sheath together during advancement.

Interference Fit is designed to ensure a snug, smooth dilator-to-sheath transition to minimize the potential for tissue trauma during vascular access.



Polyisoprene gasket is designed to be durable and enable catheter movement with leak-free hemostasis.

Integrated Suture Ring is designed to assure sheath stability.

Kink resistant strain relief – flexible design maintains sheath integrity regardless of patient movement.

Ordering Information:

REF Number	French Size	Dilator Length (cm)	Guidewire Diameter	Case Qty.
10cm length with guidewire				
<small>(Each set includes one 50cm double-straight flex-tip/J-tip guidewire and J-tip straightener)</small>				
610646	10	20*	.038	5
610650	11	20*	.038	5
610652	12	20*	.038	5
610654	13	20*	.038	5
610559	14	20*	.038	5

* Stepped Dilator

REF Number	French Size	Dilator Length (cm)	Recommended Guidewire Diameter	Case Qty.
PTCA 23cm length without guidewire				
610576	10	33*	.038	5
610577	11	33*	.038	5
610145	12	33*	.038	5
610146	13	33*	.038	5
611099	14	33*	.038	5

The Intradyn™ hemostasis introducers are designed to facilitate percutaneous introduction of catheters into the vascular system while minimizing blood loss. Refer to the Intradyn™ introducers labeling for complete instructions for use.

For more information or to place an order,
contact your B. Braun Interventional Systems Inc.
representative or call 1-877-VENA-CAV (836-2228)



Manufacturer:
B. Braun Interventional Systems Inc.
824 Twelfth Avenue
Bethlehem, PA 18018 USA
Tel: 1-877-VENA CAV (836-2228)
Fax: 610-266-3982
www.bisusa.org

Rx only
CV1093-06/11

©2011 B. Braun Interventional Systems Inc.