

PFM Medical, Inc.

Mr. Salvatore Palomares, RAC
Director of Regulatory Affairs

1815 Aston Ave, Suite 106,
Carlsbad, CA 92008, USA

pfm medical mepro gmbh

Christian Baumeister
Am Söterberg 4
66620 Nonnweiler, Germany
T +49 (0)6873 9011-42
F +49 (0)6873-9011-99
christian.baumeister@pfmmedical.com

www.pfmmedical.com

10.01.2014


To whom it may concern.

pfm medical mepro gmbh is the manufacturing facility for the Nit-Occlud® PDA device, a permanently implanted prosthesis indicated for percutaneous, transcatheter closure of small to moderate size patent ductus arteriosus with a minimum angiographic diameter less than 4 mm.

We hereby confirm, that the device is not made with natural rubber latex, i.e. it is manufactured without the use of either natural rubber latex or synthetic derivatives of natural rubber latex.

Kindest Regards,

pfm medical mepro gmbh



i. A. Christian Baumeister
Head of R&D

pfmmedical

pfm medical mepro gmbh
Am Söterberg 4 · 66620 Nonnweiler
T +49(0)6873 9011-0 · F -99
www.pfmmedical.com