

B. Braun Interventional Systems Tear Away Introducer Sheath


A Whirlwind of Options




- The B. Braun Interventional Systems Tear Away Introducer Sheath ensures a clean, even tear every time.
- Highly kink-resistant PTFE sheath material assures sheath integrity and eases device insertion.
- Ultra smooth dilator taper and dilator-to-sheath transition facilitates atraumatic device insertion.
- Unique locking system with funnel-shaped sheath entrance.

Tear Away Introducer Sheath

Tear Away Introducer Sheath Kit:

- Tear Away Introducer Sheath
- 18 Gauge x 7 cm Introducer XTW Needle*
- 10 ml Syringe
- .038" x 50 cm Double Flex Tip / 3 mm J-ST, SS Guidewire**



Ref Number	Tear Away Introducer Sheath	XTW Needle	Syringe	Stainless Steel Guidewire	Case Qty	Color Code
612600	4F x 14 cm PTFE	19Ga. x 7cm	10ml	.028" x 50cm	5	Red
612601	5F x 14 cm PTFE	18Ga. x 7cm	10ml	.038" x 50cm	5	Grey
612602	6F x 14 cm PTFE	18Ga. x 7cm	10ml	.038" x 50cm	5	Green
612603	7F x 14 cm PTFE	18Ga. x 7cm	10ml	.038" x 50cm	5	Yellow
612604	8F x 14 cm PTFE	18Ga. x 7cm	10ml	.038" x 50cm	5	Blue
612606	9F x 14 cm PTFE	18Ga. x 7cm	10ml	.038" x 50cm	5	Black
612608	10F x 14 cm PTFE	18Ga. x 7cm	10ml	.038" x 50cm	5	Pink
612609	10.5F x 14 cm PTFE	18Ga. x 7cm	10ml	.038" x 50cm	5	Light Green
612610	11F x 14 cm PTFE	18Ga. x 7cm	10ml	.038" x 50cm	5	Light Blue
612611	12F x 14 cm PTFE	18Ga. x 7cm	10ml	.038" x 50cm	5	Brown
612612	13F x 14 cm PTFE	18Ga. x 7cm	10ml	.038" x 50cm	5	Purple
612613	14F x 14 cm PTFE	18Ga. x 7cm	10ml	.038" x 50cm	5	Red
612614	15F x 14 cm PTFE	18Ga. x 7cm	10ml	.038" x 50cm	5	Grey
612615	16F x 14 cm PTFE	18Ga. x 7cm	10ml	.038" x 50cm	5	Green

*612600 has a 19 gauge x 7cm introducer XTW needle.

**612600 has a .028" x 50cm double flex tip / 3mm J-ST, SS Guidewire.

Mini Tear Away Introducer Sheath Kit:

- Tear Away Introducer Sheath
- 21 Gauge x 7 cm Introducer XTW Needle
- .018" x 45 cm Nitinol, Platinum Tip Guidewire

Ref Number	Mini Tear Away Introducer Sheath	XTW Needle	Platinum Tip Guidewire	Case Qty	Color Code
612645	4F x 7cm PTFE	21Ga. x 7cm	.018" x 45cm	10	Red
612646	5F x 7cm PTFE	21Ga. x 7cm	.018" x 45cm	10	Grey
612647	6F x 7cm PTFE	21Ga. x 7cm	.018" x 45cm	10	Green
612648	4F x 10cm PTFE	21Ga. x 7cm	.018" x 45cm	10	Red
612649	5F x 10cm PTFE	21Ga. x 7cm	.018" x 45cm	10	Grey
612650	6F x 10cm PTFE	21Ga. x 7cm	.018" x 45cm	10	Green

For more information or to place an order, contact your B. Braun Interventional Systems representative or call 1-877-VENA CAV (836-2228)

These introducers are indicated for the percutaneous introduction of diagnostic or therapeutic devices, such as catheters and pacing leads into the vasculature. Refer to the Tear Away instructions for use for relevant warnings, precautions, complications and contraindications.



LATEX
Components and packaging do not contain natural rubber latex



Manufacturer:
GALT MEDICAL CORP.
2220 Merritt Drive
Garland, TX 75041

Manufactured for:
B. Braun Interventional Systems Inc.
824 Twelfth Avenue
Bethlehem, PA 18018
Tel: 1-877-VENA CAV (836-2228)
Fax: 610-849-1334
www.bisusa.org

Rx only
CV 1076 12/12
©2012 B. Braun Interventional Systems Inc.

Mini Tear Away Introducer Sheath

Ref Number	Mini Tear Away Introducer Sheath	Mini Tear Away Introducer Sheath Length	Case Qty	Color Code
612794	4F	7cm	10	Red
612795	5F	7cm	10	Grey
612796	6F	7cm	10	Green
612797	4F	10cm	10	Red
612798	5F	10cm	10	Grey
612799	6F	10cm	10	Green
612825	7F	7cm	10	Yellow
612826	7F	10cm	10	Yellow
612830	8F	7cm	10	Blue
612827	8F	10cm	10	Blue
612831	9F	7cm	10	Black
612828	9F	10cm	10	Black
612832	10F	7cm	10	Pink
612829	10F	10cm	10	Pink

14cm Tear Away Introducer Sheath

Ref Number	Tear Away Introducer Sheath	Tear Away Introducer Sheath Length	Case Qty	Color Code
612780	4F	14cm	10	Red
612781	5F	14cm	10	Grey
612782	6F	14cm	10	Green
612783	7F	14cm	10	Yellow
612616	8F	14cm	10	Blue
612617	9F	14cm	10	Black
612619	10F	14cm	10	Pink
612620	10.5F	14cm	10	Light Green
612621	11F	14cm	10	Light Blue
612622	12F	14cm	10	Brown
612623	13F	14cm	10	Purple
612624	14F	14cm	10	Red
612625	15F	14cm	10	Grey
612626	16F	14cm	10	Green

Long Tear Away Introducer Sheath

Ref Number	Tear Away Introducer Sheath	Tear Away Introducer Sheath Length	Case Qty	Color Code
612790	4F	30cm	5	Red
612791	5F	30cm	5	Grey
612792	6F	30cm	5	Green
612793	7F	30cm	5	Yellow
612627	8F	30cm	5	Blue
612628	9F	30cm	5	Black
612629	10F	30cm	5	Pink
612630	11F	30cm	5	Light Blue
612631	12F	30cm	5	Brown
612632	13F	30cm	5	Purple
612633	14F	30cm	5	Red
612634	15F	30cm	5	Grey
612635	16F	30cm	5	Green
612850	4F	44cm	5	Red
612851	5F	44cm	5	Grey
612852	6F	44cm	5	Green
612853	7F	44cm	5	Yellow
612636	8F	44cm	5	Blue
612637	9F	44cm	5	Black
612638	10F	44cm	5	Pink
612639	11F	44cm	5	Light Blue
612640	12F	44cm	5	Brown
612641	13F	44cm	5	Purple
612642	14F	44cm	5	Red
612643	15F	44cm	5	Grey
612644	16F	44cm	5	Green