



**B. Braun Interventional Systems, Inc.**

824 Twelfth Avenue  
Bethlehem, PA 18018 USA  
Tel: 1-877-VENA CAV (836-2228)  
Fax: 610-266-3982  
[www.bbraunusa.com](http://www.bbraunusa.com)

Dear Ms. Sparling,

Thanks for taking my call today.

Further to our conversation, please consider this to be our confirmation that the Vena Tech LGM and Vena Tech LP vena cava filters are MRI compatible at 1.5T as defined by FDA.

The Vena Tech filters are non-ferromagnetic, and patients can therefore go directly into the magnet following the filter implantation procedure, if that is deemed clinically desirable.

In addition, testing of the Vena Tech LP filter at 3T has just been completed, and the results show that the filter is compatible with 3T as well.

Please contact me if you have any further questions.

Sincerely,

Paul O'Connell  
Vice-President  
Vascular Interventional Products Group  
B. Braun Medical Inc.  
800.523.9676 ext. 4818  
cell: 847.274.0097  
pager: 800.404.7103